

GSTIN: 24DCAPP9328H1ZG

Ankur Patel  
Mo:9726923885



**UMA PHARMATECH MACHINERY**

The Feature Of Biotech

**ISO 9001:2015 CERTIFIED COMPANY**

**Quality Management System**

**CE COMPLIANCE CERTIFIED**

**Sub: Quotation of Double Decker Rotary Shaker**

**Manufacturing Unit:**

**UMA PHARMATECH MACHINERY**

14, Narayan Chamber, Jogeshwari Estate, Sureliya Road  
Ahmedabad-380026, Gujarat, India

**Phone no** : + 91 -9726923885/ 8866137364

**Email** : [umapharmatechmachinerys@gmail.com](mailto:umapharmatechmachinerys@gmail.com)

[info@umapharmatechmachinery.com](mailto:info@umapharmatechmachinery.com)

**Email** <http://www.umapharmatechmachinery.com>

<http://www.fermenterbioreactor.com>

**Whats app/We Chat /IMO** : + 91-9726923885



**Skype QR Code:** This is your unique QR code for another person to scan

**Instagram:** [ankur.patel.31586](https://www.instagram.com/ankur.patel.31586)

**Bank Detail:**

**Bank Name:** Saraspur Nagarik Co-Op. Bank Ltd.

**Branch :** Vastral

**Current Account No FOR NEFT/RTGS/ECS :** 003111101004813

GSTIN: 24DCAPP9328H1ZG

PAN NO: DCAPP9328H

UAN NO: GJ01A0091839

FACTORY:2, Gajanand Industrial Hub, Behind Metrix Plaza,  
Ramol Hathijan Road, Ahmedabad, Gujarat, India

Address: 14, Narayan Chamber, Jogeshwari Estate, Sureliya  
Road, Ahmedabad-380026, Gujarat, India

Mo:9726923885

Email: [umapharmatechmachinerys@gmail.com](mailto:umapharmatechmachinerys@gmail.com)



IFS Code : GSCBOUSNCBL

GST REGISTRATION NO : 24DCAPP9328H1ZG

TRADE NAME : UMA PHARMATECH MACHINERY

Respected Sir,

This refers to our discussion regarding your requirement of Rotary Shaker  
As required by you, we are sending herewith our offer for the same.

We now request you to kindly go through our offer and if you need any further clarification, please feel free to contact us.

Thanking you and looking forward to have our long and fruitful association.

GSTIN: 24DCAPP9328H1ZG

Ankur Patel  
Mo:9726923885



# UMA PHARMATECH MACHINERY

The Feature Of Biotech



GSTIN: 24DCAPP9328H1ZG  
PAN NO: DCAPP9328H  
UAN NO: GJ01A0091839  
FACTORY:2, Gajanand Industrial Hub, Behind Metrix Plaza,  
Ramol Hathijan Road, Ahmedabad, Gujarat, India

Address: 14, Narayan Chamber, Jogeshwari Estate, Sureliya  
Road, Ahmedabad-380026, Gujarat, India  
Mo:9726923885

Email: [umapharmatechmachinerys@gmail.com](mailto:umapharmatechmachinerys@gmail.com)



## Heavy Duty Production Model

Description	Qty	Amount(INR)
<b>Rotary Flask Shaker(Double Decker)</b> <b>Platform Size( 1250 x 1250)</b> <b>Flask Arrangement</b> <b>Lower desk Capacity : 5L -25 Nos Of Flask Holding Clamp</b> <b>Upper desk Capacity : 5L -25 Nos Of Flask Holding Clamp</b>	1 NOS	Amount Without Tax & Transportation <b>1,94,900/-</b>
<ul style="list-style-type: none"> <li>• Top Sheet Made Of Mild Steel Material For Flask Mounting</li> <li>• Base Structure Made Of Mild Steel Material Painted With Red Oxide</li> <li>• Motor: Indution Gear Motor 2HP</li> <li>• Rotate With Eccentric Mechanism <b>(5 Point Technology)</b></li> <li>• <b>Speed Control VFD Drive</b></li> <li>• Outer Cover Made Of Powder Coated Mild Steel Material</li> <li>• Fixed Eccentricity <b>Stroke 30MM</b></li> <li>• This Prize Not Include Flask, clamp and rubber</li> </ul>		

### SS FLASK HOLDING CLAMP

DESCRIPTION	PRIZE PER CLAMP
<i>SS CLAMP FOR 50 ml FLASK</i>	75
<i>SS CLAMP FOR 100 ml FLASK</i>	90
<i>SS CLAMP FOR 250 ml FLASK</i>	100
<i>SS CLAMP FOR 500 ml FLASK</i>	130
<i>SS CLAMP FOR 1000 ml LASK</i>	150
<i>SS CLAMP FOR 2000 ml FLASK</i>	170
<i>SS CLAMP FOR 5000 ml FLASK</i>	200

**Validity:** 10 Days from Date of Quote.

**Supply Period:** Shall be approx. 60 to 65 Days. The delivery time will commence from the date of receipt of



*Technically and commercially clear purchase order with advance amount. Deelivery Time May Be Differ If any Problems in Manufacturing or Absences Of Workers and Engineers*

**Installation Charges:** 1. ereaction should be done by customer As Per User Manual  
2. All boarding, transportation, logging Exp in Buyer account

**PRICES:** The prices are **Ex-Factory** only.

**TAXES:** GST 18% Extra.

**PACKING:** Packing charges 3% for the instrument is extra. Based on client preference, wooden packing, plywood packing will be done at additional charges. Fumigation Charges For Export Country Is Extra Based on Service.

**FREIGHT OR TRANSPORTATION:** Extra @ actual or arrange by buyerExtra @ actual or arrange by buyer  
Machine will be dispatched on "To Pay" basis, where in the client needs to deposit the Transportation Charges at the Transporter's local office while taking the delivery. Transit Insurance is Client's responsibility. Cost would be charged extra against Handling & Forwarding Charges up to Transporter's Place at Ahmedabad.

**INSPECTION & TRAINING:** Final inspection would be carried out at our place before dispatch. We will inform the inspection date. You have to send your representative to inspect the proper functioning before dispatch. you Can testing Batch Before Dispatch Of Machine At Our Place By Using Your Equipment and instrument. Processing Parameters and Control Of Bacteria Growth are Depends On User Processing Parameters it's not our Subject. We Will Covered Hydro Test and Material Test Before Dispatch Of Machine. FAT at our works before dispatch. It is **COMPULSORY** that your technician is trained at our works for three days prior to dispatch of the machinery. A skilled technician is necessary for the same. We shall provide all the necessary training regarding the working/maintenance of the machine, the training charges will be free but the lodging and boarding and conveyance etc of our technician will be borne by you. If your operator is not trained we are not responsible for the satisfactory working of the machine. Calibration Of Ph and Do Electrode Must Have To Done By Customer Before Every Batch Start In Buyer Lab With The Use Of Suitable Buffer Solution. Accuracy Of PH Sensor  $\pm 0.5\%$   
**REJECTION CLAUSE:** If you're good self / representative fails to arrive for Final Inspection, then it will be considered that the offered machine perfectly fits in with your needs and is devoid of any manufacturing regards ever in future. Machine / goods once sold would not be taken back.

**DEMURRAGE CHARGES:** Any delay from the client's end in depositing Balance Payment after machine readiness information is passed onto him, then in such scenario we will charge 12% demurrage charges per month on the balance amount.

**VALIDITY:** The given quote is valid for 100 days only. Price change may be there beyond the validity period, subject to the market price of raw materials, availability and labour.

All Disputes Subject to Ahmedabad Jurisdiction Only

**Installation & Commissioning:** At extra cost. We will depute our service engineer(s) to commission the Instrument on your request. The charges will be confirmed on further hearing from you and paid by You. However, all Civil, Electrical and Plumbing jobs associated with the commissioning will be out of our scope. Necessary Oil/Grease and ancillaries for the commissioning to be provided by the customer. The cost of installation and commissioning will be Rs. 5,000/- per person, per day plus to and fro charges for the visit. Before the visit of the technician, please ensure that all civil and ancillary works are completed. Samples of bottles, caps, label rolls should be available in sufficient quantity. In case the above is not ready during the visit of our technician, they will go back, and RE-VISIT charges will be to your account. Customer has to provide the required samples (bottles, caps, Ampoule, Label etc) required in the trial of the machine at his cost. Quantity of the samples will be mentioned prior to the trial. Please make sure the trail samples are arranged in advance so this may not delay the delivery time. Or the delivery time scheduled in the quotation will not be valid

**Service Charges:**

Any service visit after successful completion of installation and commissioning will be charged

GSTIN: 24DCAPP9328H1ZG

Ankur Patel  
Mo:9726923885



**UMA PHARMATECH MACHINERY**

The Feature Of Biotech

@ Rs.5,000/- per person, per day, plus to and fro charges. We advise that prior to calling our engineer; your trained technician should try and troubleshoot the problem over telephone, as mostly these are small adjustment issues. Service engineers are costly and we do not want to waste either yours or our money.

**PAYMENT TERMS:**

70% deposit required to start work with P.O. Balance 30% before dispatch of machine.  
We trust that you will find our quote satisfactory and look forward to working with you. Please contact us should you have any questions at all.

FOR, UMA PHARMATECH MACHINERY



UPM

GSTIN: 24DCAPP9328H1ZG

PAN NO: DCAPP9328H

UAN NO: GJ01A0091839

FACTORY:2, Gajanand Industrial Hub, Behind Metrix Plaza,  
Ramol Hathijan Road, Ahmedabad, Gujarat, India

Address: 14, Narayan Chamber, Jogeshwari Estate, Sureliya  
Road, Ahmedabad-380026, Gujarat, India  
Mo:9726923885

Email: [umapharmatechmachinerys@gmail.com](mailto:umapharmatechmachinerys@gmail.com)