

**ISO 9001:2015 CERTIFIED COMPANY****Quality Management System****CE COMPLIANCE CERTIFIED****Sub: Quotation Of Double Wall SS Vertical Autoclave*****Buyer Detail*****Company Name :xxx****Address: xxx****Contact Number : xxx,xxx,xxx****Email : xxx,xxx****Contact Person : xx****Your Ref From :As Per Indiamart Inquiry****Your Inquiry Date: 1-April-2023****Our Offer No: UPM-2023-048****Sr No: 48****Quotation Date : 1-April-2023****Manufacturing Unit:*****UMA PHARMATECH MACHINERY***14, Narayan Chamber, Jogeshwari Estate, Sureliya Road
Ahmedabad-380026, Gujarat, India**Phone no** : + 91 -9726923885/ 8866137364**Email** : umapharmatechmachinerys@gmail.cominfo@umapharmatechmachinery.com**Email** <http://www.umapharmatechmachinery.com><http://www.fermenterbioreactor.com>

GSTIN: 24DCAPP9328H1ZG

Ankur Patel
Mo:9726923885



UMA PHARMATECH MACHINERY

The Feature Of Biotech

Whats app/We Chat /IMO : + 91-9726923885



Skype QR Code:

Instagram: [ankur.patel.31586](#)

Bank Detail:

Bank Name: Saraspur Nagarik Co-Op. Bank Ltd.

Branch : Vastral

Current Account No FOR NEFT/RTGS/ECS : 003111101004813

IFS Code : GSCBOUSNCBL

GST REGISTRATION NO : 24DCAPP9328H1ZG

TRADE NAME : UMA PHARMATECH MACHINERY

GSTIN: 24DCAPP9328H1ZG

PAN NO: DCAPP9328H

UAN NO: GJ01A0091839

FACTORY:2, Gajanand Industrial Hub, Behind Metrix Plaza,
Ramol Hathijan Road, Ahmedabad, Gujarat, India

Address: 14, Narayan Chamber, Jogeshwari Estate, Sureliya
Road, Ahmedabad-380026, Gujarat, India

Mo:9726923885

Email: umapharmatechmachinerys@gmail.com



QUOTATION

Respected Sir,

This refers to our discussion regarding your requirement of Vertical Autoclave Sterilizer(Double Wall)
As required by you, we are sending herewith our offer for the same.

We now request you to kindly go through our offer and if you need any further clarification, please feel free to contact us.

Thanking you and looking forward to have our long and fruitful association.

Description	Qty.	Prize(INR)
Double Wall Vertical Autoclave-Both Wall SS 304		
VERTICAL AUTOCLAVE – 400mm Dia X 600mm height(40 cm Dia x 60 cm Height) (TOTAL CHAMBER SIZE) Total Chamber Volume: 75(approx.) Power: 3 kw Sheet Thickness : 16 guage <ul style="list-style-type: none"> Digital Temperature Controller With Control Panel Timer For Auto Sterile 	1 Nos	+91-9726923885
VERTICAL AUTOCLAVE – 400mm Dia X 700mm height(40 cm Dia x 70 cm Height) (TOTAL CHAMBER SIZE) Total Chamber Volume: 87(approx.) Power: 3 kw Sheet Thickness : 16 guage <ul style="list-style-type: none"> Digital Temperature Controller With Control Panel Timer For Auto Sterile 	1 Nos	+91-9726923885
VERTICAL AUTOCLAVE – 470mm Dia X 580mm height(47 cm Dia x 58 cm Height) (TOTAL CHAMBER SIZE) Total Chamber Volume:100(approx.) Power: 6 kw Sheet Thickness : 2mm <ul style="list-style-type: none"> Digital Temperature Controller With Control Panel Timer For Auto Sterile 	1 Nos	Prize +91-9726923885



UMA PHARMATECH MACHINERY

The Feature Of Biotech

VERTICAL AUTOCLAVE – 470mm Dia X 725mm height(47 cm Dia x 72.5 cm Height) (TOTAL CHAMBER SIZE) Total Chamber Volume:125(approx.) Power: 6 kw Sheet Thickness : 2mm <ul style="list-style-type: none"> Digital Temperature Controller With Control Panel Timer For Auto Sterile 	1 Nos	Prize +91-9726923885
VERTICAL AUTOCLAVE – 470mm Dia X 870mm height(47 cm Dia x 87 cm Height) (TOTAL CHAMBER SIZE) Total Chamber Volume:150(approx.) Power: 6 kw Sheet Thickness : 2mm <ul style="list-style-type: none"> Digital Temperature Controller With Control Panel Timer For Auto Sterile 	1 Nos	Prize +91-9726923885
VERTICAL AUTOCLAVE – 470mm Dia X 1000mm height(47 cm Dia x 100 cm Height) (TOTAL CHAMBER SIZE) Total Chamber Volume:175(approx.) Power: 6 kw Sheet Thickness : 2.5mm <ul style="list-style-type: none"> Digital Temperature Controller With Control Panel Timer For Auto Sterile 	1 Nos	Prize +91-9726923885
VERTICAL AUTOCLAVE– 550mm Dia X 900mm height(55 cm Dia x 90 cm Height) (TOTAL CHAMBER SIZE) Total Chamber Volume:210(approx.) Power: 9 kw Sheet Thickness : 2.5mm <ul style="list-style-type: none"> Digital Temperature Controller With Control Panel Timer For Auto Sterile 	1 Nos	Prize +91-9726923885
VERTICAL AUTOCLAVE– 625mm Dia X 900mm height(62.5 cm Dia x 90 cm Height) (TOTAL CHAMBER SIZE) Total Chamber Volume:275(approx.) Power: 9 kw Sheet Thickness : 3 mm <ul style="list-style-type: none"> Digital Temperature Controller With Control Panel Timer For Auto Sterile 	1 Nos	Prize +91-9726923885



Salient Features

- Common chamber for steam and water
- chamber made of stainless steel SS 304 grade
- Outer Wall made of stainless steel SS 304 grade having 0.5mm thickness
- Non Contact Part Made Colored Mild Steel
- Lid, flange & bottom sheet also made of stainless steel. All joints argon welded.
- (Knob) Type locking system for Top Entry
- silicon gasket
- Heavy duty industrial heater.
- Design Pressure: 2 kg/cm²
- Maximum Working Pressure: 1.2Kg/cm²
- Sterilization Pressure range : 15 Psi At 121°C
- Temperature Range : 121°C
- Accessories: pressure gauge , safety spring loaded ,pressure Release valve.
- Complete unit is Hydraulically tested at (1.8 kg/cm²) pressure.
- Digital Temperature Controller With Control Panel
- Timer For Auto Sterile
- Basket Charge Extra(As Per Size)

Optional Accessories

• buzzer System	1 Nos	3600/-
• Automatic Pressure cut off	1 Nos	7900/-
• Automatic water cut off	1 Nos	7400/-

Validity: 10 Days from Date of Quote.

Supply Period: Shall be approx. 30 to 35 Days. The delivery time will commence from the date of receipt of Technically and commercially clear purchase order with advance amount. Delevery Time May Be Differ If any Problems in Manufacturing or Absences Of Workers and Engineers

Installation Charges: 1. ereaction should be done by customer
2. All boarding, transportation, logging Exp in Buyer account

PRICES: The prices are Ex-Factory only.

TAXES: GST 18% Extra.

PACKING: Packing charges 3% for the instrument is extra. Based on client preference, wooden packing, plywood packing will be done at additional charges. Fumigation Charges For Export Country Is Extra Based on Service.

FREIGHT OR TRANSPORTATION: Extra @ actual or arrange by buyerExtra @ actual or arrange by buyer
Machine will be dispatched on "To Pay" basis, where in the client needs to deposit the Transportation Charges at the Transporter's local office Before Dispatch. Transit Insurance is Client's responsibility. Cost would be charged extra against Handling & Forwarding Charges up to Transporter's Place at Ahmedabad.



UMA PHARMATECH MACHINERY

The Feature Of Biotech

INSPECTION & TRAINING: (As Per USER Manual Done By Customer) Final inspection would be carried out at our place before dispatch. We will inform the inspection date. You have to send your representative to inspect the proper functioning before dispatch. you Can testing Batch Before Dispatch Of Machine At Our Place By Using Your Equipment and instrument. We Will Covered Hydro Test and Material Test Before Dispatch Of Machine. FAT at our works before dispatch. It is **COMPULSORY** that your technician is trained at our works for three days prior to dispatch of the machinery. A skilled technician is necessary for the same. We shall provide all the necessary training regarding the working/maintenance of the machine at our site. If your operator is not trained we are not responsible for the satisfactory working of the machine.

REJECTION CLAUSE: If you're good self / representative fails to arrive for Final Inspection, then it will be considered that the offered machine perfectly fits in with your needs and is devoid of any manufacturing regards ever in future. Machine / goods once sold would not be taken back.

DEMURRAGE CHARGES: Any delay from the client's end in depositing Balance Payment after machine readiness information is passed onto him, then in such scenario we will charge 12% demurrage charges per month on the balance amount.

VALIDITY: The given quote is valid for 100 days only. Price change may be there beyond the validity period, subject to the market price of raw materials, availability and labour.
All Disputes Subject to Ahmedabad Jurisdiction Only

Service Charges:

Any service visit after successful completion of installation and commissioning will be charged @ Rs.5,000/- per person, per day, plus to and fro charges. We advise that prior to calling our engineer; your trained technician should try and troubleshoot the problem over telephone, as mostly these are small adjustment issues. Service engineers are costly and we do not want to waste either yours or our money.

PAYMENT TERMS:

70% deposit required to start work with P.O. Balance 30% before dispatch of machine.
We trust that you will find our quote satisfactory and look forward to working with you. Please contact us should you have any questions at all.

FOR, UMA PHARMATECH MACHINERY





List Of Product We Are Manufacturing

Laboratory Fermenter(1Liter to 20Liter)(Insitu or Autoclavable)

Pilot scale fermenter(50Liter to 1000L)

Production Scale Fermenter(1000Liter to 35000Liter)

Vertical Autoclave

Steam Generator(5kg/hr to 75kg/hr)(20L to 150L) steam output

Heavy Duty Shaker Machine

Incubator Shaker Machine

Water Bath Shaker

Double Cone Blender

Ribbon Blender

V Type Blender

Octagonal Blender

SS Storage Vessel

SS Mixing Vessel

SS Jacketed Processing Vessel

SS Jacketed Steam Paste Kettle

SS Filling Vessel

Dairy Storage Tank

GSTIN: 24DCAPP9328H1ZG

Ankur Patel
Mo:9726923885



UMA PHARMATECH MACHINERY

The Feature Of Biotech

Product Gallery



Laboratory fermenter



vertical Autoclave



Electric Boiler



Pilot Scale Fermenter



Production Scale SS Fermenter



Pilot Scale Fermenter



**QUALITY IS
OUR COMMITMENT**
**BUY FROM THE MANUFACTURE
SAVE YOUR MONEY**
UMA PHARMATECH MACHINERY
M: 9726923885
EMAIL: umapharmatechmachinerys@gmail.com
Address: 14, Narayan Chambers, Jogeshwari Est,
Sureliya Road, Ahmedabad-380026, Gujarat, India

Biofertilizer Plant



Fermenter



Rotary Shaker

GSTIN: 24DCAPP9328H1ZG
PAN NO: DCAPP9328H
UAN NO: GJ01A0091839
FACTORY:2, Gajanand Industrial Hub, Behind Metrix Plaza,
Ramol Hathijan Road, Ahmedabad, Gujarat, India

Address: 14, Narayan Chamber, Jogeshwari Estate, Sureliya
Road, Ahmedabad-380026, Gujarat, India
Mo:9726923885
Email: umapharmatechmachinerys@gmail.com