

# Standard Model of rotary shaker/ Orbital Shaker

## Heavy Duty Industrial Model

### Platform Size(mm) 1250 x 1250

Flask Capacity(ml)	500	1000	2000	5000
Number Of Flask	81	64	36	25

### Platform Size(mm) 1000 x 1000

Flask Capacity(ml)	500	1000	2000	5000
Number Of Flask	49	36	25	9

## Laboratory scale Model

### Platform Size(mm) 850 x 850

Flask Capacity(ml)	250	500	1000	2000	5000
Number Of Flask	49	36	25	16	9

### Platform Size(mm) 600 x 600

Flask Capacity(ml)	250	500	1000	2000	5000
Number Of Flask	25	16	9	9	4

**UMA PHARMATECH MACHINERY**

**The Feature Of Biotech**



M: +91-9726923885

Email: [info@ssfermenterbioreactor.com](mailto:info@ssfermenterbioreactor.com)

Web: [www.umapharmatechmachinery.com](http://www.umapharmatechmachinery.com)

Address: Plot No. 61, 62, Yamuna Industrial Park, Dhamatvan, Ahmedabad-382435, Gujarat, India